

## SINGLE USE DISPOSABLE COTTON ROLLS

### CONFIGURATION OVERVIEW

<b>Material:</b>	100% Cotton
<b>Characteristics:</b>	Non-Sterile, soft texture, pliable and absorbent
<b>Intended use:</b>	To adsorb oral liquid or use as wound dressing
<b>Length:</b>	38 mm $\pm$ 1.0 mm



### INSTRUCTION FOR USE

- Sterilize before use
- Choose right size according to function purpose
- Place it to mouth for absorbing liquid or use as wound dressing when necessary

### STORAGE

5 years, store it in a cool and dry place

### PRODUCT VARIATIONS

Art-No.	CORBA/CORAA26-AKZ-WH-C10	CORBA/CORAA26-AKZ-WH-C09	CORBA/CORAA26-AKZ-WH-C08
<b>Size</b>	<b>1</b>	<b>2</b>	<b>3</b>
<b>Width (mm) <math>\pm</math> 0.5 mm</b>	8	10	12
<b>PCS per bundle <math>\pm</math> 1 pcs</b>	55	40	30

### PACKAGING

CORBA26-AKZ-WH-C10/C09/C08 in boxes

CORAA26-AKZ-WH-C10/C09/C08 in bags

<b>Bundle/carton:</b>	20 x 300 g (15bundle/bag,1bag/box, 20box/ct)
<b>Carton size:</b>	41,5 x 41,5 x 35 cm
<b>Carton weight:</b>	6,5 kg
<b>Carton/pallet:</b>	24
<b>Carton/layer:</b>	4
<b>Layer/pallet:</b>	6
<b>Pallet size:</b>	80 x 120 x 230 cm
<b>Pallet weight:</b>	171 kg

<b>Bundle/carton:</b>	20X300 g (15 bundle/bag, 20bags/carton)
<b>Carton size:</b>	39 X 32 X 40 cm
<b>Carton weight:</b>	6,00 Kg
<b>Carton/pallet:</b>	30
<b>Carton/layer:</b>	6
<b>Layer/pallet:</b>	5
<b>Pallet size:</b>	80 x 120 X 220 cm
<b>Pallet weight:</b>	195 Kg